


PHILIPS

Water tank

Blue

GC092



For the storage of water in your iron

Check compatibility data

An extra-large detachable 1.4-litre tank, for the storage of water you need for ironing with steam. Allows you to iron for up to two hours without refilling!

Ease of using Philips spare parts

- Renew your product with original Philips parts

Highlights

Easy replacement

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality and 1:1 exchange.

Specifications

Replaceable part

Fits product types: GC8320/02, GC8320/07

