

PHILIPS

Steam iron

2500 series

Dripstop

Double anti-calc

GC2520



Simple, fast and effective

Longer-lasting steam performance

For great results, you want powerful steam and an iron that never gives up. With its 90 g steam boost, 2000 W for constant high steam and Double Active Calc system, this practical iron gives you value for money that lasts!

Easily removes creases

- Steam boost up to 90 g for the most stubborn creases
- Continuous steam up to 30 g/min for better crease removal
- Vertical steam for crease removal in hanging fabrics

Powerful performance

- 2000 Watt enables constant high steam output

Easy gliding on all fabrics

- Non-stick soleplate coating

Comfortable ironing

- Less refilling with extra-large 300 ml water tank

Longer lifetime

- Double Active Calc System prevents scale build-up

Highlights

2000 Watt

2000 W
max

2000 Watt enables constant high steam output.

Steam boost up to 90 g



The iron's 90 g steam boost enables you to easily remove even the most stubborn creases.

Double Active Calc System



Double Active Calc system of your Philips steam iron prevents scale build-up through anti scale pills and an easy to use calc clean function.

Specifications

Comfortable ironing

Additional comfort: 360 degree cord freedom
Anti-scale management: Double Active Calc System

Cord length: 2 m

Safe in use: Exceeds international drop test standards

Suitable for tap water

Crease removal

Continuous steam: Up to 30 gr/min

Reaching tricky areas: Pointed tip Spray

Steam Boost: Up to 90 gr/min

Vertical Steam

Easy to use

Control: Temperature ready light, Variable steam settings

Fast heat-up

Reaching tricky areas: Button groove

Smooth gliding

Soleplate: Non-stick soleplate

Technical specifications

Frequency: 50-60

Power: 2000

Voltage: 220 - 240

Weight and dimensions

Product dimensions: 29.5 x 11.6 x 19.5

Product weight: 1,2

