



Philips Azur Performer Plus
Refurbished Steam iron

Refurbished product

Steam 50 g/min; 210 g steam boost
T-ionicGlide soleplate
Safety Auto off + Anti-calc



GC4526

Faster*, Easier and Smarter

The Philips steam iron Azur Performer Plus combines powerful performance with ease of use. The built-in calc container delivers effective calc removal for long-lasting steam performance, auto steam control and T-ionicGlide soleplate .

Comfortable ironing

- Auto Steam Control for the right steam for each garment
- Innovative CordGuide
- Optimal weight for easier manoeuvrability over your garments
- Electronic shut-off when iron is left unattended
- Triple precision tip for optimal control and visibility

Calc solution

- Built-In calc container collects calc during ironing
- Self-Clean for effective calc removal

Easy gliding on all fabrics

- T-ionicGlide: our best 5-star rated soleplate

Fast and powerful performance

- Steam boost up to 210g
- 2600 W for quick heat-up and powerful performance
- Steam output of up to 50 g/min for better crease removal

PHILIPS

Highlights

Auto Steam Control



Auto Steam Control takes away the worry of selecting the amount of steam you need. Just select the temperature for the garments you're ironing and off you go.

Built-In calc container



Our specially designed and built-in calc container collects calc particles during ironing. The Self Clean process flushes the calc out of the iron to maintain consistent ironing performance day after day.

T-ionicGlide soleplate



T-ionicGlide is our best gliding and scratch-resistant, 5-star rated soleplate with integrated Titanium-Oxide layer.

Steam boost up to 210g



Steam boost up to 210g for fast removal of even the most stubborn creases.

2600 W for quick heat-up



2600 W for quick heat-up and powerful performance.

Self-Clean



In addition to the built-in calc container, the self-clean function of your soleplate will remove the most stubborn calc. A clean soleplate will improve your steam performance and ironing results.

CordGuide



The innovative CordGuide simply clicks onto your ironing board and guides the cord away while you iron.

Steam up to 50 g/min



Continuous steam output of up to 50 g/min gives you the perfect amount of steam to efficiently remove all creases.

Optimal weight



The iron has an optimal weight for moving easily over any garment, making ironing less of a hassle and making it easy to move the iron constantly between the board and its heel rest.

Refurbished Steam iron

Refurbished product Steam 50 g/min; 210 g steam boost, T-ionicGlide soleplate, Safety Auto off + Anti-calc

GC4526/87R1

Specifications

Easy to use

- Water tank capacity: 300 ml
- Power cord length: 2.5 m
- Safety auto off
- Soft grip

Scale management

- Descaling and cleaning: Built-in calc container, Self-clean

Accessories included

- CordGuide

Size and weight

- Packaging dimensions (W x H x D): 34,1 x 30 x 36 cm
- Product dimensions (W x H x D): 16 x 32,6 x 14 cm
- Total weight with packaging: 1.52 kg
- Weight of iron + base: 1.75 kg

Fast and powerful crease removal

- Continuous steam output: 50 g/min
- Power: 2600 W
- Soleplate: T-ionicGlide
- Steam boost: 210 g



Issue date 2020-10-21

Version: 1.0.1

12 NC: 0000 041 04335
EAN: 87 10103 79930 6

© 2020 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com

* Compared with an average 2200W Philips steam iron