

PHILIPS

Steam iron

PerfectCare Azur

Steam 50g/min; 190g steam boost

T-ionicGlide soleplate

Safety Auto off + Anti-calc

Safe for all ironable garments



GC4914/20



Carefree ironing, no setting required

Ironing perfected

Iron any ironable garment from silk, to linen, to cotton, to cashmere; in any order without adjusting the temperature. PerfectCare Azur delivers perfect results without the risk of burn or shine. Ironing is now simpler than ever.

One perfect setting for all your clothes

- OptimalTemp technology: perfect combination of temperature
- 100% safe on all ironable garments
- 100% fast on all fabrics, no other steam iron is faster
- 100% easy to use, no adjustment required

Comfortable ironing

- T-ionicGlide: our best 5-star rated soleplate
- Suitable for tap water with Double Active Calc Clean system
- Iron switches off automatically when left unattended
- Quattro Precision Tip for an easy reach in tricky areas
- Iron with smart light feedback indicator

Fast and powerful ironing

- 2600 W iron for fast heat-up and powerful performance
- Steam output of up to 50 g/min for better crease removal
- Steam boost up to 190 g

Highlights

OptimalTemp technology



The latest revolution in ironing to deliver the perfect combination of steam and temperature. It's made to ensure you with speedy ironing, great results on tough creases, no setting required and safe on all ironable fabrics. The perfect combination of steam and temperature because: 1) Smart Control Processor setting the right temperature 2) HeatFlow technology enabling an even steam and temperature balance.

100% safe



100% safe on all fabrics even the most delicate like silk, cashmere, wool, polyester. Independent iron testing institutes have used PerfectCare on the most sensitive ironable garments and they confirmed the excellence of ironing results.

100% fast

100% fast to iron, no sorting required. Iron all your garments with more effective steam.

100% easy to use



100% easy to use, no adjustment required. You can now iron all ironable garments one after the other, without having to wait or to adjust the iron's temperature dial.

T-ionicGlide soleplate




T-ionicGlide is Philips most premium soleplate. Delivering a new standard in glide and scratch resistance for OptimalTemp irons. The patented new coating will secure excellent results. The carefully designed shape and vent holes provide an even steam distribution for easy crease removal with the steam iron.

Suitable for tap water



PerfectCare Azur steam iron is designed to be used with tap water. Double Active Calc Clean: A smart combination of Calc Pills and a regular Self Clean reduce the scale built-up. For better ironing results we recommend to use demineralized water. This ensures a consistent steam flow from your iron.

Safety auto off



The iron shuts off automatically when left unattended. If left on it's heelrest the steam iron will switch off in 3 minutes, when left resting on the soleplate it will take 2 minutes before the iron shuts off.

Quattro Precision Tip



The Quattro Precision Tip of the steam iron allows you to reach into the most difficult to reach areas while ironing. Iron closely around the buttons with the button groove; having a speedy designed shape for better visibility; a Precision SteamBoost area for a focussed shot of steam and the slim tip for delicate ironing complete the Quattro Precision Tip.

Smart light feedback



Smart light feedback indicates the status of the iron at any time.

Specifications

Fast & powerful crease removal

Soleplate: T-ionicGlide
OptimalTEMP technology
Continuous steam output: Up to 50 g/min
Steam boost: 190 g
Power: 2600 W
Spray

Comfortable ironing

Cord length: 2.5 m

Easy to use

Water tank capacity: 350 ml
Drip stop
Safety auto off
Cord storage: Cord clip
Cord freedom (swivel): 180 degree cord freedom
Fast and easy filling
Power cord length: 2.5 m

Sustainability

Energy saving: 10% energy reduction

Calc management

Suitable for tap water
Calc clean solution: Double active calc clean

Technical specifications

Weight of iron: 1.75 kg

