

PHILIPS

ME Computerized
electric pressure cooker

Viva Collection

Press Touch

6 Liter



HD2139/05



More taste less time

This pressure cooker tenderizes foods in minutes

With a variety of cooking program, Philips electronic pressure cooker help to tenderize tough materials in minutes.

Smart and fast cooking

- Nutritional baking feature brings crispy taste
- Nutritional keep warm keeps rice fresh for 12 hours
- 0-59 mins adjustable period to keep pressure

Easy to operate

- Preset-timer for 24 hours

The ultimate in convenience

- Easy-to-read large screen digital display
- Durable and non-stick inner pot

Safe cooking

- Sensible valve locks up upper lid under high pressure
- Auto pressure release device to ensure safety during cooking

Highlights

Auto pressure release device

Auto pressure release device to ensure safety during cooking

Sensible valve

Sensible valve locks up upper lid under high pressure

Preset-timer for 24 hours



Preset-timer for 24 hours

Nutritional baking feature

Nutritional baking feature brings crispy taste

0-59 mins adjustable period to



0-59 mins adjustable period to keep pressure

Specifications

Accessories

Measuring cup
Spatula
Ladle
Soup scoop

Design specifications

Color of control panel: Black
Color(s): Champagne
Materials of main body: Color steel
Weight (incl. packaging): 6.731 kg

General specifications

14 cooking menus
Artificial Intelligence control
Backup memory when power interruption
Durable, extra thick inner pot ensures even result
Nutritional keep warm keeps rice fresh for 12 hrs
Reheat function for: fresh rice instantly
Easy-to-clean non-stick inner pot
Large LCD with clock and timer display
Detachable power cord for convenient storage

Technical specifications

Capacity: 6/12 Litres / cups
Cord length: 1.0 m
Frequency: 50 Hz
Voltage: 220 V
Wattage: 1000 W

Design

Color of main body: Champagne

Dimensions

Packaging dimensions (WxHxD):
440x350x305
Set dimension (WxHxD): 390.5x297.2x232.5

