

# PHILIPS

ME Computerized  
electric pressure cooker

Viva Collection

Press Touch

6 Liter

HD2139/65



## More taste less time

### Tenderize foods in minutes

With a variety of cooking program, Philips electronic pressure cooker help to tenderize tough materials in minutes.

#### Smart and fast cooking

- Nutritional baking feature brings crispy taste
- Nutritional keep warm keeps rice fresh for 12 hours
- 0-59 mins adjustable period to keep pressure

#### Easy to operate

- Preset-timer for 24 hours

#### The ultimate in convenience

- Durable and non-stick inner pot
- Stainless steel outer shell for easy cleaning
- Easy-to-read large screen digital display

#### Safe cooking

- Sensible valve locks up upper lid under high pressure
- Auto pressure release to ensure safety

# Highlights

## 0-59 mins adjustable period to



0-59 mins adjustable period to keep pressure

## Safety auto pressure release

In case of abnormal high pressure inside the cooker, steam will be released automatically from the pressure regulator valve to ensure safety.

## Nutritional baking feature

Nutritional baking feature brings crispy taste

## Preset-timer for 24 hours



Preset-timer for 24 hours

## Sensible valve

Sensible valve locks up upper lid under high pressure

# Specifications

## Accessories

Ladle  
Measuring cup  
Spatula

## Design specifications

Color of control panel: Black & champagne  
Color(s): Champagne metallic  
Materials of main body: Stainless steel  
Weight (incl. packaging): 6.731 kg

## General specifications

14 cooking menus  
Artificial Intelligence control  
Backup memory when power interruption  
Durable, extra thick inner pot ensures even result  
Nutritional keep warm keeps rice fresh for 12 hrs  
Reheat function for: fresh rice instantly  
Detachable and washable inner lid  
Easy-to-clean non-stick inner pot  
Handled inner pot  
Large LCD with clock and timer display

## Technical specifications

Capacity: 6 Litres / cups  
Cord length: 1.2 m  
Frequency: 50 Hz  
Voltage: 220 V  
Wattage: 1000 W

## Design

Color of main body: champagne

## Dimensions

Packaging dimensions (WxHxD):  
440x350x305  
Set dimension (WxHxD): 390.5x297.2x232.5

