

PHILIPS

Rice cooker

Healthy Variety

30 menus

10 cups



HD4777/00



Break the routine

Rice cooker with 30 varieties

The Philips' Healthy Variety rice cooker HD4777/00 come with unbeatable amount of cooking menus, Smartouch panel guiding you through cooking and a dedicated set of accessories, surprising your family with different healthy meals every day.

Smart and automatic cooking

- Artificial Intelligence control for fresh & nutritious meals
- Extra thick inner pot for more even heating

30 different healthy varieties

- 21 food menus maximize healthy varieties
- 9 rice menus for nutritious and tasty rice
- Nutritional Keep warm keeps rice fresh for 12 hours

Smartouch control makes cooking simple

- Step by step visual and audio guidance
- Advance panel design via touch navigation
- Color backlit LCD with text and animation

The ultimate in convenience

- Dedicated accessories for yoghurt, double boil and desserts
- Detachable and washable inner lid
- Inner pot with cool-touch handles for easy access

Highlights

AI control

Computerized heating mechanism and temperature control for optimal cooking result in different programs

Extra thick inner pot

Extra thick cast iron aluminum nature inside coating ensure high heating power is generated and evenly heated around the pot

9 rice menus

Rice menu includes plain rice, quick plain rice, small portion plain rice, sushi rice, brown rice, glutinous rice, mixed rice, claypot, pau fan.

21 food menus



It has 21 varieties of cooking menu including congee- plain congee thick, plain congee thin, green bean congee, millet congee, chicken congee , Pasta - with pasta sauce, Soup- vegetable soup, tom yum soup, long boil soup, Oatmeal, Yoghurt, Dessert- green bean soup, sesame paste, papaya fungus soup, lotus seed & longan soup; cake; Stew, Steam, Double Boil- soup, egg/milk custard.

Step by step guidance

Each cooking program will come with step by step visual and audio guidance during set up and cooking. User can simply follow the

instructions shown on the screen. It provides three language settings options for text guidance.

Color backlit LCD

LCD with color backlit indicates different cooking status clearly (set up in white; cooking/reheat in red; keep warm in green). It displays texts and animations to guide user through all cooking steps.

Dedicated accessories

It comes with dedicated accessories for unique cooking menus including yoghurt, sweet tofu, double boil for soup and milk custard/egg. The detachable design of steam tray can be separated as a plastic stand to cope with different cooking needs. All accessories are stackable for convenient storage.

Detachable inner lid

Wash the inner lid frequently to avoid odor. Simply remove the inner lid from the appliance and clean thoroughly.

Inner pot with handles



Easy carrying out the innerpot when you cook different varieties.

Specifications

General specifications

Artificial Intelligence control

9 rice menus

21 cooking menus

Reheat function for rice

Small portion of rice cooking

12 hours keep warm

Quick cook for plain rice

Step by step guidance

3 language setting options

Touch sensor control panel

Color backlit LCD display

Easy-to-program timer

Durable, extra thick inner pot ensures even result

Easy clean non-stick innerpot

Handled inner pot

Advanced lid design

Dedicated cooking accessories

Detachable inner lid

Spillover prevention vent

Swing handle for easy carrying

Detachable power cord for convenient storage

Backup memory in power failure

Water level indicator

Accessories

Rice scoop

Soup scoop

Measuring cup

Plastic steam tray

Double boil pot

Yoghurt cup

Technical specifications

Capacity: 1.8/10 Litres / cups

Inner pot capacity: 5 litres

Cord length: 1.0 m

Voltage: 220 V

Frequency: 50 Hz

Wattage: 825 (China) W

Design specifications

Materials of main body: PP

Weight appliance: 5.18 kg

Weight (incl. packaging): 7.5 kg

Color(s): White body, Black panel

Width: 287 mm

Height: 289.7 mm

Depth: 249.4 mm

