

PHILIPS

Water tank

HD5093



For the water in your espresso machine

Please check specifications for compatible products

This water container fits Saeco, Saeco and Gaggia Accademia espresso machines. Its handle makes it easy to detach for refilling and cleaning.

Ease of using Philips spare parts

- Easily renew your product with original Philips parts

Highlights

Easy replacement of product parts

From time to time your product needs a facelift, and with Philips consumer replacement parts to renew your product, it's never been so easy! All this with guaranteed Philips quality.

Specifications

Replaceable part

Fits Coffee Product: HD8953/01, HD8953/09, HD8953/11, HD8953/19, HD8953/21, HD8954/01, HD8954/09, HD9712/01, HD9712/11

