

PHILIPS

Saeco

Manual Espresso

Class



HD8325



Authentic Italian Espresso every day

With pressurized filter Crema

Philips-Saeco manual espresso machine offers enthusiasts of traditional preparation everything for a perfect daily espresso. The patented pressurized filter Crema guarantees a long lasting delicious crema, at every use.

Authentic Italian Espresso

- Pressurized filter Crema
- 15 bar pump
- Cup warmer

Simplicity and speed of use

- Grind coffee and Easy Serving Espresso (E.S.E.)
- Pannarello to froth milk as you like
- Ergonomic in day by day operations

Every day use

- Long lasting premium materials

Highlights

Pressurized filter Crema



This special crema filter guarantees a long lasting delicious crema

Grind coffee and Easy Serving



For your convenience you can choose between ground coffee and Easy Serving Espresso (E.S.E.) pads.

Pannarello



With the pannarello you can froth your milk exactly according to your preferences. Just dip the pannarello in the milk, while slowly stirring

the milk with a rotary movement. Froth the milk as long as you want for more or less foamy mikfroth.

15 bar pump



High pressure ensures that the coffee beans entire aroma is always fully extracted.

Long lasting premium materials



Philips Saeco Manual Espresso features high quality and durable housing material

Cup warmer



This feature allows you to not only store your cups and glasses right on the machine, it also pre-heats them letting the aroma unfold even better and guaranteeing a long-lasting crema

Ergonomic in day by day



Refilling coffee powder or water, emptying the filter or the drip tray. All compartments are directly accessible for maximum convenience.

Specifications

Features

Brita Filter compatible: Optionable
Cupholder
Cupwarmer
Filterholder: Pressurized "Crema"
Hot water / Steam Nozzle
Interface: Rotary Switch
Kit Pods E.S.E
Pannarello: Stainless Steel
Simultaneous brewing

Technical specifications

Boiler material: Inox

Technical data

Power: 950 W
Pump Pressure: 15 bar

Logistic data

Pallet number of pieces: 50
Pallet size: 80,0 x 120,0 x 205,0 cm
Pallet weight (EU): 235 kg

Water tank

Water tank capacity: 1 l

