

PHILIPS

Steamer

Viva Collection

9 L 900 W

Manual timer

HD9120/93



Discover fuller flavours

Add your own taste with the flavour booster

The Flavour Booster provides your food with great aromas of herbs and spices. You can prepare favourite soup, stew, porridge, rice and more in the XL steaming bowl. Just use the timer and off you go!

Discover variety of fuller flavours

- Flavour Booster adds more taste with delicious herbs & spices
- XL steaming bowl for soup, stew, rice and more

Only a little effort

- Dishwasher-safe parts
- 60-minute timer
- External water inlet

Smart features

- Removable 500 ml drip tray is easy to pour from and clean
- Cord storage
- Stackable steaming tiers
- Removable bottoms to steam larger pieces of food

Highlights

Flavour Booster



The unique Flavour Booster of the Philips steamer adds delicious aroma of herbs and spices, bringing even more taste to steaming. Simply pop your favourite herbs and spices into the booster, and let steam do the rest. The heat from the steam releases delicate aromas from the herbs and spices, which thoroughly infuse the food with their mouth-watering flavours.

XL steaming bowl



With the XL steaming bowl you can prepare more varied tasty meals for your family. XL steaming bowl for soup, stew, rice and more.

60-minute timer



60-minute timer with ready signal and auto shut off.

External water inlet



External water inlet for refilling reservoir during use.

Stackable steaming tiers



Stackable steaming tiers use less space for storing.

Dishwasher-safe parts



Dishwasher-safe parts make cleaning easy.

Cord storage

Cord storage keeps kitchen surface neat and tidy.

Removable 500 ml drip tray

Removable 500 ml drip tray is easy to pour from and clean

Removable bottoms

Need to steam a large piece of food, like a whole chicken, or any dish that's too big for a single steaming compartment? Simply unclip the bottom of each basket to create a large single space. The steamer will do the rest for you, steaming thoroughly and consistently from top to bottom.

Specifications

Design specifications

Material: Plastic
Color(s): Black with silver accents

Accessories included

XL steaming bowl: 2.5 L
Egg rack

General specifications

Non-slip feet
Power-on light
Cool-touch handgrips

Water level indicator
Overheat and dry boil protection
Number of tiers: 3 pcs

Weight and dimensions

Weight appliance: 2.3 kg
Product dimensions (W x H x D): 364 x 450 x 228 mm

Technical specifications

Power: 900 W
Voltage: 220-240 V
Frequency: 50/60 Hz
Cord length: 100 cm
Capacity (max): 2.5 / 2.6 / 3.5 L
Capacity water tank: 1.1 L

