

PHILIPS

Steamer

Avance Collection

9 L 1350 W

Digital timer, heating turbo

Aroma Infuser, meal program

Plastic/Metal, White

HD9171/00



The purest taste of steamed food

Each food cooked at right time for a tastier meal

The meal program of the Philips steamer lets you set the right cooking time for different foods. Auto timer tells you when to add each basket at the right order, no more over or under cooked food. Whole tasty meal, all cooked at right times

Discover variety of fuller flavours

- Aroma Infuser adds more taste with delicious herb & spices
- Recipe booklet

Smart features

- Stackable steaming tiers
- Water refill indicator

Only a little effort

- Dishwasher-safe parts
- Digital timer
- Keep Warm function keeps food ready for serving


Highlights

Aroma Infuser



The unique Aroma Infuser of the Philips steamer adds delicious aroma of herbs and spices, bringing even more taste to steaming. Simply pop your favourite herbs and spices into the booster, and let steam do the rest. The heat from the steam releases delicate aromas from the herbs and spices, which thoroughly infuse the food with their mouth-watering flavours.

Recipe booklet



Recipe booklet packed with expert tips and tasty treats.

Stackable steaming tiers



Stackable steaming tiers use less space for storing.

Water refill indicator



Water refill indicator signals when water level is low.

Dishwasher-safe parts



Dishwasher-safe parts make cleaning easy.

Digital timer



Digital timer with ready signal and auto shutoff.

Keep Warm function



Keep Warm function keeps your food warm till it's served.

Specifications

Design specifications

Material: Plastic and metal
Color(s): White and stainless steel

Accessories included

Egg rack
Rice/soup/food container: 2.5 L

General specifications

Non-slip feet
Power-on light
Cool-touch handgrips
Water level indicator
Number of tiers: 3 pcs
Overheat and dry boil protection
Cord storage

Technical specifications

Power: 1350 W
Voltage: 220-240 V
Frequency: 50/60 Hz
Cord length: 100 cm
Capacity (max): 2.6/ 2.8/ 3.2 L
Capacity water tank: 1.5 L

