

PHILIPS

Airfryer XXL

Accessory kit

1x pizza tray

1x recipe booklet



HD9953/00



Pizza Master Kit

Accessories and tips to master Airfryer Pizza

With this special Philips Airfryer Pizza master kit, you can become a true Pizza Master. Without any pre-heating needed, you can enjoy your favorite pizza in the Airfryer in just 8 minutes!

Mouthwatering Pizza made easy

- Bake your favorite 26cm in only 8 min. using the Pizza tray
- Get inspired by the tips, tricks from the recipe booklet

Only a little effort

- Dishwasher-safe parts for easy cleaning

Highlights

Pizza Tray

Now you can use your Airfryer to bake your favorite pizza using the Pizza Tray. It will bake 26cm home-made pizza or even frozen pre-made ones as fast as only 8 minutes.

Inspiring booklet



Booklet with chef's tips, beginner's recipe and Airfryer cooking times is included.

Easy cleaning and storing

You can safely put the pizza tray and tongs in your dishwasher making them even easier to re-use.

Specifications

Accessories included

Pizza tray
Recipe booklet

Design and finishing

Material of main body: Plated steel

Product Compatibility

Compatible with: Philips Airfryer XXL:
HD963x and HD965x

Weight and dimensions

Dimensions of product (LxWxH):

354x240x165 mm

Weight of product: 1.55 kg

Dimensions of packaging (LxWxH):

375x250x185 mm

Weight incl. packaging: 1.68 kg

Sustainability

Packaging: > 90% recycled materials

User manual: 100% recycled paper

