

Health and Well-being range 2001



HF 340

HF 350

Picture not yet available

HF 341

Picture not yet available

HF 346

Picture not yet available

HF 351

Picture not yet available

HF 380



HF 305

Picture not yet available

HF 306

Picture not yet available


HF 318



HF 319



HF 321



HF 365

VITAL AND BODY SIGNS

	HF 340 Bathroom scales 200g	HF350 Bathroom scales 100g	HF341 Bathroom scales available sept 2001	HF346 Bathroom scales available sept 2001	HF351 Bathroom scales available sept 2001	HF380 Bathroom scales available sept 2001	HF305 Upper arm blood pressure monitor	HF306 Upper arm blood pressure monitor available sept 2001	HF318 Wrist blood pressure with pouch available sept 2001	HF319 Wrist blood pressure monitor	HF321 Wrist blood pressure with pouch	HF365 Electronic digital thermometer
Features												
Display	LED	LED	LED	LCD	LCD	LCD	LCD	LCD	LCD	LCD	LCD	LCD
Power	9 V alkaline, not included	9 V alkaline, not included	9 V alkaline, not included	2xCR2032 Lithium included	2xCR2032 Lithium included	2xCR2032 Lithium included	4 x AA Alkaline (Philips Power Life) included	4 x AA Alkaline (Philips Power Life) included	2 x AA Alkaline (Philips Power Life) included	2 x AA Alkaline (Philips Power Life) included	2 x AA Alkaline (Philips Power Life) included	1.55 V, LR41 included
Resolution	200 g	100 g	200 g	200 g	100 g	100 g	1 mm Hg	1 mm Hg	1 mm Hg	1 mm Hg	1 mm Hg	0.01 °C
Features	blue top anti-slip surface large clear display large sturdy feet	white top anti-slip surface large clear display large sturdy feet	large clear display mat white top weight lock auto off	large clear display pastel green top lithium battery weight lock auto off	large clear display mat white top with silver ring lithium battery weight lock auto off	large clear display frosted glass top lithium battery weight lock auto off	Fuzzy logic controlled automatic inflation Memory Systolic, diastolic and pulse rate displayed simultaneously Storage pouch Personal registration card	Fuzzy logic controlled automatic inflation 60 memories Systolic, diastolic and pulse rate displayed silent pumping Storage pouch Personal registration card	Automatic inflation Pulse rate no memories Systolic, diastolic and pulse rate displayed simultaneously Personal registration card	Automatic inflation Pulse rate 60 memories Systolic, diastolic and pulse rate displayed simultaneously Personal registration card	Automatic inflation Pulse rate 60 memories Systolic, diastolic and pulse rate displayed simultaneously Storage pouch Personal registration card	Waterproof Buzzer Memory translucent blue protecting case
Accessoires												
Logistic data												
Dimensions single packaging (Fbox) lwxhx (mm)	332 x 68 x 366	332 x 68 x 366	332 x 68 x 366	332 x 68 x 366	332 x 68 x 366	336 x 53 x 323	212 x 92 x 184	187 x 92 x 170	110 x 93 x 155	110 x 93 x 155	110 x 93 x 155	55 x 23 x 150
EAN Fbox 8710103+	102977	103349	150701	150725	150749	150763	103363	150688	150664	137108	137122	103417
												display box for 12 pieces 174 x 100 x 153 EAN 8710103 103424
Dimensions A box packaging lwxhx (mm)	374 x 280 x 346	374 x 280 x 346	374 x 280 x 346	374 x 280 x 346	374 x 280 x 346	not available	376 x 220 x 198	376 x 195 x 182	226 x 186 x 159	226 x 186 x 159	226 x 186 x 159	527 x 378 x 334
EAN A box	102984	103356	150718	150732	150756	150770	193370	150695	150671	137115	137139	103431
Products/ A box	4	4	4	4	4	4	4	4	4	4	4	20 x 12 = 244



PHILIPS

Let's make things better.



HP 5225



HF 3304



HF 3305

RELAXATION

HP5225
Foot bath massager + whirl

Features

Display	-
Power	230 V, 70 Watt
Resolution	-
Features	Two massage speeds Keep warm setting Whirl Cord storage on bottom
Accessoires	-

Logistic data

Dimensions single packaging (Fbox) l x w x h (mm)	245 x 375 x 420
EAN Fbox	8710103+ 024774
Dimensions A box packaging l x w x h (mm)	520 x 390 x 440
EAN A box	055440
Products/ A box	2

BRIGHT LIGHT

HF 3304

HF 3305

Features

Wattage (W)	72	72
Lamps	2xPhilips Bright Light (PL-L 36 W)	2xPhilips Bright Light (PL-L 36 W)
2500 Lux level (cm)	60	60
Electronic ballast	●	●
Dimmer	-	●
Desk position	●	●
Wall position	-	-
Floorstand	-	-

Logistic data

Packaging size (cm)	36 x 21.5 x 62.2	36 x 21.5 x 62.2
Product weight (kg)	± 3.1	±3.1
Pallet quantity	30	30

International range (note: not all models will be available in all countries)