

PHILIPS

Digital Scale



HF8000



Give your body an accurate voice

Fast Response

Philips electronic scale gives you an accurate weight measurement. It handles weights up to 150kg with a precision of 100g. Stylish, convenient and easy to use.

Professionally accurate measurement

- Measures up to 150 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

Fast Response

- Measurement results within seconds
- Instantly ready to use: simply step on

Specifications

Features

Easy to clean
Immediately ready to use
Anti-slip safety

Technical Specifications

Maximum weight: 150 kg
Minimum weight: 10 kg
Resolution: 100 g
Unit of measure: kg / st / lb
Gravity compensation: 6 zones
Technology: 4-sensor
Display: LED (can also be used in dark places)

Color: white

Dimensions: 326 x 326 x 30 mm

Batteries: 4 x LR6 AA 1.5V

Accessories

Manual

