

PHILIPS

Electronic scale



HF8003



Give your body an accurate voice

Weight memory

Philips electronic scale gives you an accurate weight measurement. It handles weights up to 150 kg with a precision of 100 g. Stylish, convenient and easy to use.

Professionally accurate measurement

- Measures up to 150 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

Track weight changes

- Displays your previous weight

Easy readability

- Large and high contrast display for easy readability

Specifications

Features

Weight memory
Easy-to-read display
Easy to clean
Immediately ready to use
Anti-slip safety

Technical specifications

Maximum weight: 150 kg
Minimum weight: 10 kg
Resolution: 100 g
Unit of measure
Gravity compensation: 6 zones
Technology: 4-sensor
Display: LCD

Color(s): white

Dimensions: 326 x 326 x 30 mm

Batteries: 4 x LR6 AA 1.5V

Accessories

Batteries: included
Carpet feet: 4 pieces
Manual

