



Philips
Electronic scale



HF8003

Give your body an accurate voice

Weight memory

Philips electronic scale gives you an accurate weight measurement. It handles weights up to 150 kg with a precision of 100 g. Stylish, convenient and easy to use.

Professionally accurate measurement

- Measures up to 150 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

Track weight changes

- Displays your previous weight

Easy readability

- Large and high contrast display for easy readability

PHILIPS
sense and simplicity

Specifications

Features

- Weight memory
- Easy-to-read display
- Easy to clean
- Immediately ready to use
- Anti-slip safety

Technical specifications

- Maximum weight: 150 kg
- Minimum weight: 10 kg
- Resolution: 100 g
- Unit of measure

- Gravity compensation: 6 zones
- Technology: 4-sensor
- Display: LCD
- Color(s): white
- Dimensions: 326 x 326 x 30 mm
- Batteries: 4 x LR6 AA 1.5V

Accessories

- Batteries: included
- Carpet feet: 4 pieces
- Manual



Issue date 2012-03-09

Version: 2.1.1

12 NC: 0000 000 00000

© 2012 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com