

Philips Electronic scale



HF8003

# Give your body an accurate voice

### Weight Memory

Philips electronic scale gives you an accurate weight measurement. It handles weights up to 150kg with a precision of 100g. Stylish, convenient and easy to use.

#### Professionally accurate measurement

- Measures up to 150 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- · Adjusts for gravity differences around the world

#### Track weight changes

· Displays your previous weight

#### **Easy readability**

· Large and high contrast display for easy readability



Electronic scale

#### **Specifications**

#### **Features**

- Weight memory
- · Easy-to-read display
- Easy to clean
- · Immediately ready to use
- Anti-slip safety

## Technical Specifications • Maximum weight: 150 kg

- Minimum weight: 10 kg

- Resolution: 100 g
  Unit of measure: kg / st / lb,
  Gravity compensation: 6 zones
  Technology: 4-sensor
  Display: LCD
  Colours white

- Colour: white
- Dimensions: 326 x 326 x 30 mm
  Batteries: 4 x LR6 AA 1.5V

#### **Accessories**

- Batteries: included
- Carpet feet: 4 pieces
- Manual



Issue date 2009-02-11

Version: 2.1

© 2009 Koninklijke Philips Electronics N.V. All Rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips Electronics N.V. or their respective owners.

www.philips.com