

# PHILIPS

## Digital scale



HF8005



## Give your body an accurate voice

### Weight management

Philips digital scale gives you accurate weight measurement up to 180 kg. Its advanced weight management features make it easy to keep track of changes in your weight and help you on your way to good general health.

#### **Professionally accurate measurement**

- Measures up to 180 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

#### **Weight management**

- Shows your weight history over time
- Assists in working towards and maintaining your ideal weight
- Supports multiple users

#### **Easy readability**

- Large and high contrast display for easy readability

#### **Durable, easy care design**

- Durable and easy to clean stainless steel design

# Specifications

## Features

Weight management  
Weight history  
Multiple users  
Easy-to-read display  
Durable design  
Easy to clean  
Immediately ready to use

## Technical specifications

Maximum weight: 180 kg  
Minimum weight: 10 kg  
Resolution: 100 g  
Unit of measure  
Gravity compensation: 6 zones  
Technology: 4-sensor  
Display: LCD  
Color(s): Grey (brushed stainless steel)

Material: Stainless steel  
Dimensions: 298x298x40 mm  
Batteries: 4 x LR6 AA 1.5V

## Accessories

Batteries: included  
Carpet feet: 4 pieces  
Manual

