



Steam generator iron

- Max 5 bar pump pressure
- 1.1-l detachable water tank
- 180 g steam boost

HI5916/26

Faster ironing

With extra steam

This steam generator makes your ironing experience faster thanks to extra powerful continuous steam, and more convenient thanks to a 1.1 L detachable water tank.

Benefits

Compact and convenient

- 1.1-l detachable water tank
- Lightweight and compact for easy use and storage

Fast and powerful

- Calc-clean container included — no cartridges, no extra cost
- Powerful steam for ultimate crease removal
- Ceramic soleplate for durability and easy gliding

PHILIPS

Features

Detachable 1.1 l water tank



The water tank can be removed to be refilled easily under the tap. The large filling hole makes refilling fast. Thanks to the 1.1-l water tank capacity, you can iron for up to 2 hours continuously without any refills.

Smart descaling



Our built-in descaling system, Smart Calc Clean, reminds you when you need to descale. It also includes a container to make descaling easy. This means no cartridges are required and there are no additional costs.

Ultra-powerful steam



Strong, continuous steam tackles even the thickest fabrics with ease. Watch stubborn creases melt away with an extra steam boost just where you need it. This extra steam is perfect for vertical steaming to refresh clothes and curtains, too.

Compact and convenient



The lightweight, compact size is perfect for storage and fits conveniently on your ironing board. Exclusive ProVelocity technology makes our steam generators smaller and more compact than ever.

Ceramic soleplate



Our durable ceramic soleplate glides well on all garments. It is non-stick, scratch-resistant and simple to keep clean.

Specifications

Please note that this is a pre-sales leaflet. The contents of this leaflet reflect the best of our knowledge per date and country mentioned above. The contents of this leaflet are subject to change without notice. Philips does not accept any liability as to the contents of this leaflet.

Fast crease removal

Continuous steam	up to 100 g/min
Power	2400 W
Pressure	max 5 bar pump pressure
Steam boost	180 g
Vertical steam	Yes
Voltage	220-240 V

Easy to use

Soleplate name	Ceramic
Water tank capacity	1100 ml
Detachable water tank	Yes
Hose length	1.6 m
Power cord length	1.6 m
Refill anytime during use	Yes
Tap water suitable	Yes
Extra-large filling hole	Yes

Scale management

Descaling and cleaning	Smart Calc Clean
Descaling reminder	Yes no cartridges needed, no extra

Accessories included

Calc clean container	Yes
----------------------	-----

Size and weight

Product dimensions (W x H x D)	36.2 x 27 x 26.3 cm
Weight of iron	1.1 kg
Weight of iron + base	3.2 kg

Guarantee

2 year worldwide guarantee	Yes
----------------------------	-----

Technical specifications

Heating time	2 min.
--------------	--------

Packaging dimensions

Height	28.00 cm
Width	39.60 cm
Depth	23.60 cm
Net weight	2.70 kg

Gross weight 3.67 kg
EAN 08710103864790
Number of products included 1
Country of origin CN

Harmonised Systems Code 851640



data subject to change
2018, October 26
Version: 4.1.1
EAN: 08710103864790

© 2018 Koninklijke Philips N.V.
All rights reserved.

Specifications are subject to change without notice. Trademarks are the property of Koninklijke Philips N.V. or their respective owners.

www.philips.com