

PHILIPS

InfraPhil infrared
lamp

150 W



HP3616



Relieves muscular pains

Deep, focused infrared warmth

Deep, focused warmth from infrared lamp

Relieves muscle pain

- Deep penetrating Infrared warmth

Extra-focused treatment

- 150 Watt infrared lamp with extra focus

Easy to use

- Adjustable angle
- On/off switch

Easy to store

- Cord storage

Highlights

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

150 Watt infrared lamp with extra focus



150W

The Philips Infrared lamp has been optimised for pain relief in areas covering 20 x 30 cm for example shoulder, elbow, calf or neck. The concentric rings help to focus the beams on a specific area.

Adjustable angle



The infrared lamp can easily be positioned up to 40 degrees backwards.

On/off switch



Easy on/off switch at the back of the infrared lamp, so the appliance does not have to be unplugged after each treatment.

Cord storage



The infrared lamp has a cord storage holder at the back of the appliance. Simply wrap the cord around it and it can easily be stored again.

Specifications

Easy to use

Adjustability: 0–40 degrees (backwards)
On/off switch

Easy storage

Cord storage

Technical specifications

Power: 150 W

Voltage: 220/230 or 240/250 V V

Frequency: 50 Hz

Cord length: 180 cm

Insulation: Class II (double isolation)

Type of lamps: PAR 38E, 150 W + prismatic rings for more focus

Lifetime of lamp: 750 sessions of 10 minutes

Weight and dimensions

Product dimensions: 22.0 x 13.0 x 18.5 (H x W x D) cm

Product weight: 1 kg

F-box dimensions: 22.5 x 14.0 x 20.4 (H x W x D) cm

F-box weight: 1 kg

No. of F-boxes in A-box: 6 pcs

A-box dimensions: 24.0 x 28.8 x 62.4 (H x W x D) cm

A-box weight: 6.4 kg

Qty. on Euro pallet: 192 pcs

Logistic data

CTV code: 884361601000

Country of origin: Germany

Safety

IEC certified: Complies to IEC 60335

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

