

PHILIPS

InfraCare infrared
lamp

650W

Half body treatment

63.5cm high



HP3641



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared warmth is very comfortable and penetrates deep into the skin, stimulates blood circulation and gives effective relief for half-body areas of 60x40 cm.

Effective pain relief

- Deep penetrating Infrared warmth

Half-body treatment

- 650 Watt infrared halogen lamp

Easy to use

- Adjustable angle with rotation

Treatment control

- Digital timer with automatic shut-off

Personalized warmth comfort

- Dimmer with 3 warmth settings

Highlights

650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle with rotation



The InfraCare HP 3643 is easy to use for half-body-treatment and for different situations (e.g. lying on a bed or couch) due its maximum flexibility. Next to the height adjustment, the lamphouse can easily be positioned to

optimally treat the right body part, by moving the lamphouse into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Digital timer with automatic shut-off



Treatment time can easily be controlled by setting the advised or preferred number of minutes. The digital timer automatically counts down and automatically shuts off at the end of the session.

Dimmer with 3 warmth settings



The infrared warmth offered by the InfraCare can easily be adjusted to ones personal preference. Three different levels of warmth intensity can be chosen, making the treatment as comfortable as possible.

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Effective pain relief

Deep penetrating warmth

Half-body treatment

Treatment area: 60x40 cm

Infrared halogen lamp: 650 W

Treatment control

Dimmer with 3 warmth settings

Digital timer: with auto shut-off

Pause function

Easy to use

2 handgrips: for easy transport and positioning

Extendable height: 45-80 cm

Digital user interface: for feedback on time and warmth intensity

On/off switch

Easy positionable lamphouse (left/right): 90

Easy positionable lamphouse (up/down): 45

Technical specifications

Power: 650 W

Voltage: 230 V

Frequency: 50/60 Hz

Lifetime of lamps: 500 hour(s)

Cord length: 300 cm

Insulation: Class II (double isolation)

Weight and dimensions

Product dimensions: 29.0x63.5x29.0

(WxHxD) cm

Product weight: 5.5 kg

F-box dimensions (=A-box): 33.9x80.0x36.5

(WxHxD) cm

F-box weight (=A-box): 6.5 kg

Qnt. on Euro pallet: 18 pcs

Logistic data

CTV code: 884364101000

Country of origin: Germany

Safety

IEC certified: Complies to IEC 60601 and 60335

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

