

PHILIPS

InfraCare infrared lamp

650W

Half body treatment

124cm high



HP3643



Effective pain relief

Deep penetrating infrared warmth

Philips InfraCare gives effective pain relief for muscles and joints. Its 650 W infrared lamp is very comfortable and the infrared warmth penetrates deep into the skin, stimulates blood circulation and gives relief for half-body areas of 60x40 cm.

Effective pain relief

- Deep penetrating Infrared warmth

Half-body treatment

- 650 Watt infrared halogen lamp

Easy to use

- Adjustable angle with rotation
- Easily reaches neck and shoulders

Stimulated blood circulation

- Stimulates blood circulation



asimpleswitch.com

Highlights

650 Watt infrared halogen lamp



The Philips InfraCare HP 3643 has been optimised for effective relief of pains in areas covering 60 x 40 cm, such as the entire back, both shoulders and neck or thigh. The InfraCare features innovative infrared halogen lamp technology. The special optics, filter and the powerful halogen lamp of 650 Watt have been optimised, enabling half-body treatment (60 x 40cm), and much more comfortable warmth, as the infrared warmth is more equally divided over the treatment area.

Adjustable angle with rotation



The InfraCare HP 3643 is easy to use for half-body-treatment and for different situations (e.g. lying on a bed or couch) due its maximum flexibility. Next to the height adjustment, the lamphouse can easily be positioned to optimally treat the right body part, by moving the lamphouse into a horizontal or vertical position (left/right 90 degrees) and by tilting it up or down (45 degrees).

Easily reaches neck and shoulders



The InfraCare HP3643 can easily reach neck shoulders to give optimal warmth treatment of sore muscles. The height of the lamp house can be extended from 60 to 120 cm, so it offers a convenient position whether you are lying on a bed or sitting on a chair.

Stimulates blood circulation

Infrared light and warmth penetrates deep into your skin, stimulating blood circulation and dilating blood vessels. Furthermore, infrared light speeds up the metabolic processes and the removal of waste substances from the body.

Deep penetrating Infrared warmth



Infrared light has been proven to effectively relieve aches and pains caused by muscular pains and stiff joints. The comfortable warmth created by infrared light penetrates deep into your skin, stimulates blood circulation and warms your muscles. As your muscles are soothed, they automatically loosen up and relax. Because it makes the tissue more pliable, warmth also reduces stiffness and makes joints more flexible.

Specifications

Effective pain relief

Deep penetrating warmth

Half-body treatment

Treatment area: 60x40 cm

Infrared halogen lamp: 650 W

Easy to use

On/off switch

Adjustability: 0-40 degrees (backwards)

2 handgrips: for easy transport and positioning

Extendable height: 62-127 cm

Easy positionable lamphouse (left/right): 90

Easy positionable lamphouse (up/down): 45

Technical specifications

Power: 650 W

Frequency: 50-60 Hz

Lifetime of lamps: 500 hour(s)

Cord length: 300 cm

Insulation: Class II (double insulation)

Voltage: 220-230 V, Taiwan 110 V

Weight and dimensions

Product dimensions: 29x124x29 cm

Product weight: 7.2 kg

F-box dimensions: 34.5x60x39 cm

F-box weight: 9.1 kg

F-box dimensions (=A-box): 34.5x60x39 cm

F-box weight (=A-box): 9.1 kg

Qnt. on Euro pallet: 24 pcs

Logistic data

CTV code: 884364301000

Country of origin: Germany

Safety

IEC certified: Complies to IEC 60601 and 60335

15 minutes auto shut-off

Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC



asimpleswitch.com

Philips Green Logo

Philips Green Products can reduce costs, energy consumption and CO2 emissions. How? They offer a significant environmental improvement in one or more of the Philips Green Focal Areas – Energy efficiency, Packaging, Hazardous substances, Weight, Recycling and disposal and Lifetime reliability.