

# PHILIPS

## Epilator

### Satinelle Ice

HP6492/11



## New ceramic epilator from Philips Norelco

Ice and massage for a more gentle hair removal

Thanks to the soothing properties of ice and an integrated massage system, the new Satinelle Ice significantly reduces the negative sensation normally associated with epilation.

#### Advanced pain reduction system

- Integrated skin cooler
- Pivoting cooling head
- Integrated massaging system

#### Unique ceramic epilating head

- Hypo-allergenic ceramic epilating head
- Curved epilating head
- Higher rotation speed

#### Additional benefits

- For smooth shaving of underarms and bikini line
- Exfoliating body puff
- New ergonomic design
- Two speed settings

# Highlights

## Integrated skin cooler

Cools and soothes the skin during and immediately after epilation, making hair removal more gentle than ever.

## Pivoting cooling head

Follows the body's contours during epilation and immediately calms any negative sensation on the skin.

## Integrated massaging system

Works with the skin cooler to distract from the sensation of epilation.

## Hypo-allergenic ceramic epilating head

Better skin care by reducing the risk of allergic skin irritation.

## Curved epilating head

For improved grip on hair, ensuring that up to 95% of hairs are removed at the root in one stroke.

## Higher rotation speed

For more effective hair removal.

## Exfoliating body puff for less ingrown hairs

Reduces ingrown hair.

## New ergonomic design

For easier handling.

## 2 speed settings

Speed 1 for fine hair, speed 2 for thicker hair.

# Specifications

## Features

Ceramic epilation system  
Pivoting ice cooler  
Active massage system  
2 speed settings

## Technical specifications

Power source: AC (Mains)  
Voltage: 100–240V (50–60 Hz) Automatica  
Voltage Selection V  
Motor: DC-Motor 14V  
Power consumption: 6 W  
Number of discs: 25  
Number of catching points: 24  
Pulling actions/second speed 1: 880  
Pulling actions/second speed 2: 1080  
RPM speed 1: 2200 min  
RPM speed 2: 2700 min

## Accessories

Detachable shaving head  
Insulation sleeve  
Exfoliation body puff  
Luxury storage pouch  
Cleaning brush

## Weight and dimensions

Packaging design: Blister  
F-box dimensions: 6.2"x4.3"x10.4" (width x depth x height) mm  
F-box weight: 1.1 lb  
A-box dimensions: 10.4"x12.5"x7.3" (width x depth x height) mm  
A-box weight: 3.3 g  
No. of pieces per A-box: 2 pcs

