

# PHILIPS

Nutri wonder  
straightener

StraightCare  
Prestige

ActiveRelease technology

Natural oils + Vitamin E

Premium styling performance

Hair becomes softer, shinier  
and stays strong\*



HP8383/00



## Beautifully straight. Feels wonderfully nourished. Reinventing straightening with infusion covers

Straighten and protect your hair in one go, leaving it feeling wonderfully nourished. Nutri wonder comes with replaceable plate covers, proven to release nutrient-rich ingredients- Argan oil, Camellia oil, Avocado oil and Vitamin E – onto your hair as you straighten.

### Care and protect

- ActiveRelease technology – scientifically proven to release natural ingredients
- The power of natural ingredients
- Heat protection for every hair strand
- Replace the covers every month for the best results
- Ceramic coating for smooth gliding and damage prevention

### Beautifully styled hair

- Extra care and control with 3 temperature settings
- Up to 200°C for perfect results
- You can choose to style your hair with or without covers

### Ease of use

- Instant heat-up in just 15 sec
- Easy cover-release buttons on the side of the straightener
- Swivel cord to prevent tangled wires
- Safe storage of the product with soft heat-protection cap

# Highlights

## ActiveRelease technology



The plate covers contain a lightweight blend of Argan, Camellia and Avocado oils, combined with Vitamin E, which using the straightener's heat, are gradually released onto the hair as you straighten, leaving you with a cashmere-soft and radiant finish\*\*.

## Natural ingredients



Nutri wonder brings together a lightweight blend of ingredients used for centuries by women around the globe. The formula has been carefully crafted to suit all hair types. Camellia oil is a great source of Omega 3, 8 and 9. It's traditionally used for protection and miniaturisation of hair. Argan oil is a powerful mix of Vitamins E and F and other vital nutrients. It's well known for its antioxidant and nourishing properties. Avocado oil is rich in Vitamins A, D and E and proteins. It is famous for its nourishing and moisturising effect on hair.

## Heat protection



Nutri wonder straightens with the utmost care. The oil-releasing plate covers shield your hair from direct heat whilst the smooth surfaces

minimise friction and protect the health of your hair. When the plate covers are placed correctly, the 'Care' sign will appear on the display to signal that the straightener is ready to use.

## Replace the covers every month



Each set of plate covers lasts one month with average use, which is 10 straightening sessions. When the delicate flowery aroma and decorative design on top of the plates fades, it's time for new ones.

## Variable Temperature Control



Three temperature settings from 150°-200° allow you to adjust the styling temperature to suit your hair. Choose low for thin hair, medium for normal hair and high for thick hair.

## Up to 200°C



High temperature enables you to change the shape of your hair and get the perfect style you want.

## 15 sec heat-up time




No more long waits for your straightener to be ready. Nutri wonder heats up in just 15 seconds.

## Straighten with/without cover

You can choose to style your hair with or without covers, thanks to floating ceramic plates.

## Ceramic plates



The smooth ceramic-coated plates prevent damage to your hair during styling with smooth gliding.

## Easy cover release



Covers are easy to insert and detach thanks to release sliding buttons on the side of the straightener.

# Specifications

## Accessories

2 FREE plate covers included  
Soft heat protection cap

## Features

Heat/Speed settings: 3  
Coating of heated parts: Ceramic  
Auto shut-off: after 60 min

Swivel cord  
Storage hook  
Ready-to-use indicator

## Caring technologies

ActiveRelease technology  
Natural oils + Vitamin E

## Technical specifications

Variable temperature: 150–200°C  
Number of heat settings: 3  
Cord length: 1.8 m  
Heater type: High-performance heater  
Type of temperature control: Digital LCD  
Plate size: 22.5 x 90 mm

## Service

2-year worldwide guarantee

