

# PHILIPS


6 Nutri wonder plate covers+2 as a gift

**StraightCare  
Prestige**

ActiveRelease technology

Natural oils + Vitamin E

Premium styling performance



HP8394/00



## Beautifully straight. Feels wonderfully nourished. Reinventing straightening with infusion covers

Straighten and protect your hair in one go, leaving it feeling wonderfully nourished. Nutri wonder replaceable plate covers are proven to release nutrient-rich ingredients- Argan oil, Camellia oil, Avocado oil and Vitamin E – onto your hair as you straighten.

### Care and protect

- ActiveRelease technology – scientifically proven to release natural ingredients
- The power of natural ingredients
- Heat protection for every hair strand
- Replace the covers every month for the best results

### Beautifully styled hair


- You can choose to style your hair with or without covers

### Ease of use

- Easy cover-release buttons on the side of the straightener

# Highlights

## ActiveRelease technology



The plate covers contain a lightweight blend of Argan, Camellia and Avocado oils, combined with Vitamin E, which using the straighteners' heat, are gradually released onto the hair as you straighten, leaving you with a cashmere-soft and radiant finish\*\*.

## Natural ingredients



Nutri wonder brings together a lightweight blend of ingredients used for centuries by women around the globe. The formula has been carefully crafted to suit all hair types. Camellia oil is a great source of Omega 3, 8 and 9. It's traditionally used for protection and miniaturisation of hair. Argan oil is a powerful mix of Vitamins E and F and other vital

nutrients. It's well known for its antioxidant and nourishing properties. Avocado oil is rich in Vitamins A, D and E and proteins. It is famous for its nourishing and moisturising effect on hair.

## Heat protection



Nutri wonder straightens with the utmost care. The oil-releasing plate covers shield your hair from direct heat whilst the smooth surfaces minimise friction and protect the health of your hair. When the plate covers are placed correctly, the 'Care' sign will appear on the display to signal that the straightener is ready to use.

## Replace the covers every month



Each set of plate covers lasts one month with average use, which is 10 straightening sessions. When the delicate flowery aroma and decorative design on top of the plates fades, it's time for new ones.

## Straighten with/without cover

You can choose to style your hair with or without covers, thanks to floating ceramic plates.

## Easy cover release



Covers are easy to insert and detach thanks to release sliding buttons on the side of the straightener.

# Specifications

## Caring technologies

ActiveRelease technology

Natural oils + Vitamin E

