



Philips
shaving head cleaning fluid

HQ100/01

For a clean shave

Thoroughly cleans the shaving heads

Thoroughly cleans the shaving heads

- Excellent hygienic performance
- Replaceable filter

Philips shaving head cleaner

- Not suitable for 6600, 7100, 7200, 8100, 9100 series
- Suitable for both 3-head & 2-head Philips shavers
- Including 150 ml cleaning fluid

PHILIPS

Specifications

Highlights

Suitable for 600, 6000, 7000, 8000 series only

Not suitable for 6600, 7100, 7200, 8100, 9100 series

Excellent hygienic performance

Excellent hygienic performance

Replaceable filter

Replaceable filter

Suitable for both 3-head & 2-head Philips shavers

Suitable for both 3-head & 2-head Philips shavers

Including 150 ml cleaning fluid

Including 150 ml cleaning fluid



Issue date 2015-01-27

Version: 4.1.1

12 NC: 0000 000 00000

© 2015 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com