

# PHILIPS

## jet Clean cleaning solution

Cleans and lubricates

Cool Breeze scent

HQ200/50



## Keep a clean shave

### For maximum performance

A clean shaver ensures maximum shaving performance. Use your Philips Jet Clean solution HQ200/50 with your Jet Cleaning system for a thorough cleaning of your Philips shaving heads.

#### **Keeps your shaver in perfect condition**

- Ensures maximum shaving performance

#### **Cleans and lubricates your shaving heads**

- For a thorough clean

#### **New Cool Breeze Scent**

- Leaves a fresh scent on your shaving heads

# Highlights

## Jet Clean solution

Jet Clean solution cleans and lubricates your shaving heads.

## Cleans thoroughly

After cleaning, the shaving heads of your Philips electric shaver smell, look and perform as new.

## New Cool Breeze Scent

Fresh scent that will remain on your shaving head, so you can enjoy the perfect fresh and clean shave of your electric shaver.

# Specifications

## Content

**Capacity:** 300 ml

**Content sufficient for:** Up to ~1 month of use

**Packaging:** 1 bottle

## Cleaning

**Usage:** Change the shaving cleaning liquid every month

