

PHILIPS

jet Clean cleaning solution

Cleans and lubricates

Cool Breeze scent

HQ200/50



Keep a clean shave

For maximum performance

A clean shaver ensures maximum shaving performance. Use your Philips Jet Clean solution HQ200/50 with your Jet Cleaning system for a thorough cleaning of your Philips shaving heads.

Keeps your shaver in perfect condition

- Ensures maximum shaving performance

Cleans and lubricates your shaving heads

- For a thorough clean

New Cool Breeze Scent

- Leaves a fresh scent on your shaving heads

Highlights

Jet Clean solution

Jet Clean solution cleans and lubricates your shaving heads.

Cleans thoroughly

After cleaning, the shaving heads of your Philips electric shaver smell, look and perform as new.

New Cool Breeze Scent

Fresh scent that will remain on your shaving head, so you can enjoy the perfect fresh and clean shave of your electric shaver.

Specifications

Content

Capacity: 300 ml

Content sufficient for: Up to ~1 month of use

Packaging: 1 bottle

Cleaning

Usage: Change the shaving cleaning liquid every month

