

PHILIPS

jet Clean cleaning solution

Cleans and lubricates

Cool Breeze scent

Value pack

HQ203/51



Keep a clean shave

For maximum performance

A clean shaver ensures maximum shaving performance. Use your Jet Clean solution with your Jet Clean system for a thorough cleaning of your shaving heads

Value Pack

- 3 pack Jet Clean Solution

Keeps your shaver in perfect condition

- Ensures maximum shaving performance

Cleans and lubricates your shaving heads

- For a thorough clean

New Cool Breeze Scent

- Leaves a fresh scent on your shaving heads

Highlights

Value Pack

3 pack Jet Clean Solution, to enjoy longer a fresh and clean shave.

Jet Clean solution

Jet Clean solution cleans and lubricates your shaving heads.

Cleans thoroughly

After cleaning, the shaving heads of your Philips electric shaver smell, look and perform as new.

New Cool Breeze Scent

Fresh scent that will remain on your shaving head, so you can enjoy the perfect fresh and clean shave of your electric shaver.

Specifications

Cleaning

Usage: Change the liquid after each cleaning procedure

Content

Capacity: 900 ml

Content sufficient for: Up to ~ 3 months of use

Packaging: 3 bottles

