

Philishave

4000 series MICRO ACTION

HQ 4830



▼123456799▲

Philishave



is better.



PHILIPS

This Micro Action shaver with a rechargeable and mains power supply features a 2 star luxury with an attractive Forest Green finishing

shaving system





Shaving system features connected to Micro Action

- Microgroove for a closer shave
- Lift & Cut system lifts hair before cutting it off
- 9 personal comfort settings for the sensitive and less sensitive skin
- long omnidirectional slots catch and shave long and short stubble
- three 'Micro Action' shaving heads (HQ 4)

Shaving system features connected to the Unique Philishave Rotary system

- 45 razor sharp blades for an extra fast shave
- Silent Power for a silent yet powerful motor
- unique rotating system for a vibration free operation
- hair collecting chamber allows hygenical collections of hairs shaven off





power system



rechargeable + mains

- rechargeable for cordless shaving
- mains operation possible for shaving when the battery is empty
- 1 hour charging time allows quick charging
- up to 30 minutes cordless shaving time for 10 days cordless shaving
- quick charge in 6 minutes for 1 cordless shave
- worldwide automatic voltage selection between 100-240V AC 50-60 Hz
- green nickel-hydride battery contains no cadmium, mercury and lead



luxe




electronic featuring

- 1 LED charging light indicates connected to the mains (green)
- on/off switch with integrated lock

cosmetics and accessories

- Forest Green finishing
- pop-up trimmer for moustache and sideburns
- pouch for storage and protection during travelling
- coiled cord, cleaning brush and shaving head protection cap



packaging

fancy box

- weight shaver 212 grams
- dimensions fancy box 134 x 66 x 190 mm
- weight fancy box 430 grams
- dimensions master carton 275 x 200 x 216 mm
- weight master carton 3000 grams
- master carton contains 6 pieces
- EAN code 87 10103 075899
- CTV code 88 54 830 16000

blister

- weight shaver 212 grams
- dimensions blister 170 x 90 x 265 mm
- weight blister 472 grams
- dimensions master carton 468 x 296 x 186 mm
- weight master carton 3152 grams
- master carton contains 6 pieces
- EAN code 87 10103 075486
- CTV code 88 54 830 33000