

PHILISHAVE

400 series MICRO ACTION


HQ 484



rechargeable
and mains



luxe



PHILIPS

PHILISHAVE

400 series MICRO ACTION

HQ 484

this rechargeable shaver combines an innovative, daring design with extrovert colours and the excellent shaving performance of the Micro Action shaving system.


shaving system




two 'Micro Action' shaving heads, model HQ 4
Micro Groove brings the cutters even closer to the skin
Silent Power silent yet powerful motor




individual floating heads follow the contours of the face
Lift & Cut system lifts each hair before cutting it off
 adjustable shaving heads allow optimum setting for every skin sensitivity
 long omnidirectional slots catch long and short stubble
 self sharpening blades for years of optimum performance
 unique rotating system for vibration free operation




power system




rechargeable for cordless shaving
 mains operation possible for shaving when the battery is empty
 1 hour charging time allows quick recharging
 30 minutes cordless shaving time for 1-2 weeks cordless shaving
 world-wide automatic voltage selection between 100-240V AC 50-60 Hz
 'green' nickel-hydrate battery contains less than 0.01% cadmium, mercury and lead
 built-in charger eliminates the need for a separate charging unit



luxury features



pop-up trimmer for moustache and sideburns
 charging light indicates that the shaver is connected to the mains
 on/off switch with integrated lock prevents accidental operation in transit
 deluxe travel pouch with cord compartment for storage and protection during travelling
 coiled cord, cleaning brush and shaving heads protection cap



packaging

weight shaver 158 grams
 fancy box dimensions 117 x 72 x 169 mm, weight 375 grams
 master carton dimensions 238 x 219 x 184 mm, weight 2470 grams, contains 6 pcs
 EAN code: 87 10103 04361 4
 CTV code: 8851 484 16000

