


Philishave

5000 series REFLEX ACTION

HQ 5465



Philishave



...was better.



PHILIPS

This Reflex Action shaver with a mains only power supply features a 4 star luxury with a sophisticated Seaweed Metallic Green finishing

shaving system



Shaving system features connected to Reflex Action

contour following shaving system for a closer and more comfortable shave

Lift & Cut system lifts hair before cutting it off
9 personal comfort settings for the sensitive and less sensitive skin
long omnidirectional slots catch and shave long and short stubble
three 'Reflex Action' shaving heads (HQ 5)

Shaving system features connected to the Unique Philishave Rotary system

45 razor sharp blades for an extra fast shave
Silent Power for a silent yet powerful motor
unique rotating system for a vibration free operation
hair collecting chamber allows hygenical collections of hairs shaven off



power system




mains

mains operating shaver
worldwide automatic voltage selection between 100-240V AC 50-60 Hz



luxe



electronic featuring

electronic push button switch with full metal stainless steel surrounding

cosmetics and accessories

Seaweed Metallic Green finishing
pop-up trimmer for moustache and sideburns
luxurious compact cassette with cord compartment for storage and protection during travelling
coiled cord, cleaning brush and shaving head protection cap



packaging

fancy box

weight shaver 150 grams
dimensions fancy box 133 x 87 x 190 mm
weight fancy box 468 grams
dimensions master carton 283 x 200 x 277 mm
weight master carton 3093grams
master carton contains 6 pieces
EAN code 87 10103 074724
CTV code 88 55 465 01000

blister

weight shaver 150 grams
dimensions blister 170 x 90 x 265 mm
weight blister 450 grams
dimensions master carton 564 x 184 x 281 mm
weight master carton 3300 grams
master carton contains 6 pieces
EAN code 87 10103 075295
CTV code 88 55 465 33000