

Philishave

5000 series REFLEX ACTION

HQ 5830



Philishave



Shaves better.



PHILIPS

This Reflex Action shaver with a rechargeable and mains power supply features a 2 star luxury with an attractive Aubergine Black finishing

shaving system




Shaving system features connected to Reflex Action

contour following shaving system for a closer and more comfortable shave

Lift & Cut system lifts hair before cutting it off
9 personal comfort settings for the sensitive and less sensitive skin
long omnidirectional slots catch and shave long and short stubble
three 'Reflex Action' shaving heads (HQ 5)

Shaving system features connected to the Unique Philishave Rotary system

45 razor sharp blades for an extra fast shave
Silent Power for a silent yet powerful motor
unique rotating system for a vibration free operation
hair collecting chamber allows hygienical collections of hairs shaven off






power system



rechargeable + mains

rechargeable for cordless shaving
mains operation possible for shaving when the batteries are empty
1 hour charging time allows quick charging
up to 70 minutes cordless shaving time for 3 weeks cordless shaving
quick charge in 6 minutes for 1 cordless shave
worldwide automatic voltage selection between 100-240V AC 50-60 Hz
low voltage boat or car charging allows maximum travel convenience
green nickel-hydride batteries contains no cadmium, mercury and lead



luxe



electronic featuring

2 LED battery charging (green) and battery reserve control (red)
soft touch slide switch with integrated lock

cosmetics and accessories

Aubergine Black finishing
pop-up trimmer for moustache and sideburns
pouch with cord compartment for storage and protection during travelling
coiled cord, cleaning brush and shaving head protection cap



packaging

fancy box

weight shaver 204 grams
dimensions fancy box 133 x 71 x 190 mm
weight fancy box 437 grams
dimensions master carton 278 x 199 x 230 mm
weight master carton 2907 grams
master carton contains 6 pieces
EAN code 87 10103 074694
CTV code 88 55 830 16000

blister

weight shaver 204 grams
dimensions blister 170 x 90 x 265 mm
weight blister 504 grams
dimensions master carton 564 x 184 x 281 mm
weight master carton 3624 grams
master carton contains 6 pieces
EAN code 87 10103 075271
CTV code 88 55 830 33000