

# Philishave

5000 series REFLEX ACTION

HQ 5853



Philishave



as better.



# PHILIPS

# Philishave

5000 SERIES REFLEX ACTION

# HQ 5853

This Reflex Action shaver with a rechargeable and mains power supply features a 3 star luxury with a sportive Spicy Orange finishing

## shaving system




Shaving system features connected to Reflex Action

contour following shaving system for a closer and more comfortable shave

Lift & Cut system lifts hair before cutting it off  
9 personal comfort settings for the sensitive and less sensitive skin  
long omnidirectional slots catch and shave long and short stubble  
three 'Reflex Action' shaving heads (HQ 5)

Shaving system features connected to the Unique Philishave Rotary system

45 razor sharp blades for an extra fast shave  
Silent Power for a silent yet powerful motor  
unique rotating system for a vibration free operation  
hair collecting chamber allows hygienical collections of hairs shaven off





## power system




### rechargeable + mains

rechargeable for cordless shaving  
mains operation possible for shaving when the batteries are empty  
1 hour charging time allows quick charging  
up to 70 minutes cordless shaving time for 3 weeks cordless shaving  
quick charge in 6 minutes for 1 cordless shave  
worldwide automatic voltage selection between 100-240V AC 50-60 Hz  
low voltage boat or car charging allows maximum travel convenience  
green nickel-hydride batteries contains no cadmium, mercury and lead



## luxe



### electronic featuring

5 LED battery charging control (green)  
LED battery reserve control (red)  
soft touch slide switch with integrated lock

### cosmetics and accessories

Spicy Orange finishing with sportive detail colouring  
pop-up trimmer for moustache and sideburns  
luxurious pouch with cord compartment for storage and protection during travelling  
coiled cord, cleaning brush and shaving head protection cap



## packaging

### fancy box

weight shaver 210 grams  
dimensions fancy box 133 x 75 x 190 mm  
weight fancy box 498 grams  
dimensions master carton 275 x 199 x 246 mm  
weight master carton 3273 grams  
master carton contains 6 pieces  
EAN code 87 10103 074632  
CTV code 88 55 853 16000

### blister

weight shaver 210 grams  
dimensions blister 170 x 90 x 265 mm  
weight blister 510 grams  
dimensions master carton 564 x 184 x 281 mm  
weight master carton 3660 grams  
master carton contains 6 pieces  
EAN code 87 10103 075233  
CTV code 88 55 853 33000