

Philishave

Coolskin
3 headers
2 headers



Philishave range 2002



PHILIPS

Let's make things better.



COOL SKIN

	HQ 6757	HQ 6756/6755	HQ 6735	HQ 6715	HQ 673	HQ 671
Moisturising Shaving Emulsion						
Optimised for use with shaving emulsion	●	●	●	●	●	●
Dry use possible	●	●	●	●	●	●
Shaving emulsion cartridges included	5	5	5	5	5	5
Free gel samples	5	5	3	1	1	1
Washable with water	●	●	●	●	●	●
Shaving system						
3D contour following	●	●	●	●	●	●
Super Lift & Cut	●	●	●	●	●	●
Individual Floating heads	●	●	●	●	●	●
Corrosion free (CCC) shaving heads	●	●	●	●	●	●
Silent power	●	●	●	●	●	●
Long, omnidirectional slots	●	●	●	●	●	●
45 Razor sharp blades	●	●	●	●	●	●
Shaving unit release knob	●	●	●	●	●	●
Detachable flip top bracket with fixation	●	●	●	●	●	●
Floes Transparent	●	●	●	●	●	●
Ornamental plate plastic coloured	●	●	●	●	●	●
Shaving system name on deco caps	●	●	●	●	●	●
Vibration free rotating system	●	●	●	●	●	●
Power system						
Rechargeable mains	-	-	-	-	-	-
Rechargeable (deductive charging)	●	●	●	●	●	●
Power plug	-	-	-	-	●	●
Charging stand	●	●	●	●	-	-
High voltage circuitry in power plug and low voltage in shaver and stand	●	●	●	●	●	●
Charging time	1 hr	1 hr	8 hr	8 hrs	8 hr	8 hrs
Total shaving time	40'	40'	40'	40'	30'	30'
Automatic Voltage selection	●	●	●	●	●	●
12-24 Volts charging	-	-	-	-	-	-
Number of batteries required	-	-	-	-	-	-
Luxury level						
Colours	Exclusive	Luxe	Regular	Basic	Regular	Basic
Lacquer	Innovative finish	Perfect grip	Grip metallic	-	Perfect grip	-
Charging indicator	●	●	●	●	●	●
Battery low indicator	●	●	●	●	●	●
Battery full indicator	●	●	-	-	-	-
Click-on trimmer unit, optimised for moustache and side-burns	●	●	●	-	●	-
Hanging cord	●	●	●	●	-	-
Deluxe pouch with storage possibility for 2 additive cartridges	-	-	-	●	-	-
Basic travel pouch	-	-	-	-	●	●
Travel cassette	-	-	●	-	-	-
Luxe cassette	Limited Edition	●	-	-	-	-

International range (note: not all models will be available in all countries)