

PHILIPS

Pasta Maker Cookie
Accessory

Avance Collection

Cookie kit

HR2455/05



Enjoy fresh and healthy cookies with no effort

Cookie kit for the Philips noodle maker

Want to enjoy fresh cookies yet have healthy concerns? Thanks to Philips pasta maker cookie kit, you can make low fat healthy cookie at home with minimum effort now.

Low fat healthy cookies

- Special recipes to make low fat healthy cookies

Great convenience

- Mixing, kneading and extruding in one go

Highlights

Low fat healthy cookies

The kit comes with special recipes to make healthy cookies with 50% less butter and sugar than the regular recipe.

Fully automatic process

After choosing the program and pressing the start button, the machine automatically mixes, kneads and extrudes the cookies, and it only takes 10 minutes!

Specifications

Accessories

Included: Cleaning tool

