

PHILIPS

Food processor

350 W

2.1 l bowl



HR7605/10



The versatile space saver

Easy to store food processor

This Philips Food Processor is a space-saving, time-saving multifunctional answer to all of your kitchen problems. With this one appliance you can handle over 15 different jobs in a fraction of the time it would normally take.

Powerful precision

- Powerful 350 Watt motor for the toughest ingredients
- Easily performs 15+ functions

Compact

- All accessories fit in the bowl for compact storage
- Small enough to fit easily in any cupboard

Highlights

Ultra compact



Small enough to fit easily in any cupboard

MicroStore



Microstore allows you to store all your accessories inside the bowl for easy and compact storage.

15+ functions



The Food Processor comes with 5 different accessories that enable it to perform 15+ different functions. The dishwasher-safe accessories are: a kneading tool for mixing batter and for heavy kneading; a stainless steel chopping knife for preparing meat and vegetables; metal disc inserts for handling medium shredding, medium slicing and granulating; and an emulsifying disc for preparing food like whipped cream and mayonnaise.

Specifications

Accessories

Number of inserts/discs: 3

Shredding insert medium: For medium shredding of vegetables and fruit

Slicing insert medium: For medium slicing of vegetables and fruits

Granulating insert medium: For granulating of potatoes, cheese

Emulsifying disc: Whipping, Whisking, Emulsifying application

Stainless steel chopping knife: chopping, blending, kneading, pureeing, mixing

Design specifications

Housing: PP

Switch and tool holder: POM

Bowl, cover, pusher: SAN

Emulsifying disc: ABS

Metal knife and inserts: stainless steel and ABS

Technical specifications

Motor: Universal, radio/TV interference suppressed

Power: 350 W

Cake batter capacity: 250 g

Bowl content: 2.1 l

Voltage: 220-240 V

Frequency: 50-60 Hz

Colour(s): White with mineral green accents

Safety: Automatic resettable cut-off protecting the motor

Speeds: 1+ pulse

Cord length: 75 cm

Country of origin

Hungary

