



Philips NIVEA
Shaving conditioner

HS800

Moisturizing shaving conditioner

Only for Philips NIVEA FOR MEN HS8000 series

The formula is enriched with Chamomile and Vitamins to condition your skin and support a healthy look.

Moisturizing and soothing shave

- NIVEA FOR MEN conditioner dispensing system

Easy shaving conditioner refill

- Integrated cartridge refill system

PHILIPS

NIVEA
FOR MEN

Specifications

Accessories

- Shaving conditioner refill can: for Nivea for Men

Shaver 8000 series

Highlights

Integrated cartridge refill system

Simply pump the NIVEA FOR MEN shaving conditioner from the refill can directly into the Philips NIVEA FOR MEN shaver.

Conditioner dispensing system

NIVEA FOR MEN moisturizing shaving conditioner is dispensed directly onto your skin