

Software history (HTS5563, 5583, 5593)

This firmware is only applicable for devices with the following serial number: **PT1Axxxxxxxxxx**

If you device serial number is **PT2Axxxxxxxxxx**, please proceed to go HTS5563M2, HTS5583M2 or HTS5593M2 product page for device firmware upgrade

The latest software solves the issues below:

NOTE: Each release contains all solutions identified in the earlier version.

Version 1.59.00

- Fixed network error during CDDA playback
- Improved compatibility with Smartphone (MyRemote) app
- Some fixes for NetTV services
- Some models unable to wake up at correct source mode

Version 1.52.00

- Improved Audio Sync
- Fixed system bugs

Version 1.47.10

- Improved system performance