

# PHILIPS

## sonicare

Standard sonic  
toothbrush heads

ProResults


3-pack

Standard size

Click-on

All-around cleaning

HX6013



## Superior performance.\* Guaranteed quality. Superior performance at a superior value

As one of our flagship toothbrush heads, Philips Sonicare ProResults is perfect for new or legacy Sonicare users who simply desire that authentic Sonicare cleaning experience at an unbelievable value.

### **Thorough clean feeling in your mouth**

- Removes up to 2 times more plaque than a manual toothbrush

### **Designed to optimize performance**

- Optimized Philips Sonicare performance
- Reminder bristles ensure your most effective clean
- Engineered to maximize sonic motion

### **Proven to improve oral health**

- Part of a better oral healthcare routine

### **A brush head that fits multiple handles**

- Click-on design for simple brush head placement

# Specifications

## Items included

**Brush heads:** 3 ProResults standard

## Compatibility

**Brush head system:** Click-on

**Suitable for these models:** HealthyWhite+, 2 Series plaque control, 3 Series gum health,

DiamondClean, EasyClean, FlexCare, FlexCare Platinum, FlexCare+, for Kids, HealthyWhite

## Design and finishing

**Bristle stiffness feel:** Regular

**Color:** White

**Reminder bristles:** Blue bristle color fade away

**Size:** Standard

## Health benefits

**Plaque removal:** Removes up to 2x more plaque\*

## Quality and performance

**Replacement:** Every 3 months

**Tested:** for optimal usage

