

PHILIPS

sonicare

Standard sonic
toothbrush heads

C1 ProResults

3-pack

Standard size

Click-on

HX6013/31



Superior performance.* Guaranteed quality. Superior performance at a superior value

As one of our flagship toothbrush heads, the Philips Sonicare C1 ProResults is perfect for new or existing Philips Sonicare users who want our authentic cleaning experience at an unbelievable value.

For exceptional care of teeth and gums

- Up to 3x more plaque removal than a manual toothbrush

An easy addition to your oral care routine

- Get the most effective plaque removal, every day

Innovative technology

- Philips Sonicare's advanced sonic technology

Designed to optimize performance

- Tested to meet your oral health needs
- Click-on design for simple brush head placement

Highlights

Say goodbye to plaque



Densely-packed, high-quality bristles remove up to 3x more plaque than a manual toothbrush. And the contoured bristles fit the shape of your teeth, so you can easily clean those hard-to-reach areas.

Reminder bristles



Brush heads become less effective after 3 months of use, but we'll remind you before this happens. Just monitor your blue replacement bristles, and when they turn white you'll know it's time for a fresh brush head.

Easy to click-on



Your C1 ProResults brush head fits perfectly with any Philips Sonicare toothbrush handle, except Philips One and Essence. Simply click on and off for easy replacement and cleaning.

The only true sonic technology



Philips Sonicare's advanced sonic technology pulses water between teeth, and its brush strokes break up plaque and sweep it away for an exceptional daily clean.

Brushing with peace of mind



All Philips Sonicare brush heads are safe and gentle on teeth and gums. During tests, each brush head is put through its paces to ensure durability and an exceptional performance every time you brush.

Specifications

Compatibility

Brush head system: Click-on
Suitable for these models: 2 Series plaque control, 2 Series plaque defense, 3 Series gum health, DiamondClean, DiamondClean Smart, EasyClean, Essence+, FlexCare, FlexCare Platinum, FlexCare Platinum Connected, FlexCare+, for Kids, HealthyWhite, HealthyWhite+, ProtectiveClean

Design and finishing

Bristle stiffness feel: Medium
Color: White
Reminder bristles: Blue bristle color fade away
Size: Standard

Items included

Brush heads: 3 C1 ProResults standard

Quality and performance

Replacement: Every 3 months
Tested: for optimal usage

Health benefits

Plaque removal: Removes up to 2x more plaque*

