

PHILIPS

sonicare

Tongue cleaning
starter kit

TongueCare+

1 tongue brush

1 tongue spray



HX8071/23



The only tongue cleaner with Sonicare strength

The only tongue cleaner with Sonicare strength

The new Philips Sonicare TongueCare+ tongue brush and spray transform a Sonicare toothbrush into the most advanced breath management solution. Kills oral malodor for all-day fresh breath** Reduces oral malodor up to 78% better***

Innovative technology


- Designed specifically for the tongue
- A powerful combination

An easy addition to your oral care routine

- A convenient, healthy habit


Highlights

Design



In healthy patients, up to 80% of oral malodor comes from the tongue. But most patients don't clean it, and those who do typically use a manual toothbrush, which isn't designed to clean the tongue.

Dual-action technology



Philips Sonicare TongueCare+ pairs sonic technology and the antibacterial strength of Philips Sonicare BreathRx Tongue Spray to

target the primary cause of oral malodor.

Innovative tongue brush with 240 elastomer MicroBristles to clean deep around papillae. 31,000 vibrations per minute help break up tongue coating and sweep away debris and bacteria. Alcohol-free tongue spray with cetylpyridinium chloride (CPC) and Zytex™ to neutralize VSC gasses, in a thick, concentrated formula that lubricates the tongue for a gentle clean.

Convenience and compliance



The TongueCare+ tongue brush attaches to most Sonicare power toothbrush handles in just one click, making it a simple addition to your patients' oral care routine. The compact shape helps patients comfortably access difficult-to-reach areas on the back of the tongue, where most oral malodor bacteria live.

Specifications

Items included

Antibacterial tongue spray: 2oz
Tongue brush: 1

Ingredients

Antibacterial tongue spray: 0.09% CPC, Zytex

Expiration

After opened: 6 months
Shelf life: 36 months

