

Philips Sonicare
Rechargeable sonic
toothbrush



HX9553/02

Custom Care system

4 ways of customised cleaning in one product

2 brush head sizes and 2 cleaning speeds

- Standard for a complete clean, compact for precision clean
- Thorough everyday clean, massages and stimulates gums

Customisable cleaning

- Quadpacer 30-second interval timer encourages even brushing
- Easy-start to acclimate to the Sonicare experience

Superior cleaning technology

- Dynamic fluid action reaches deep between teeth
- Removes plaque in hard-to-reach areas
- Clinically proven safe and gentle
- Healthier gums

PHILIPS
sonicare
the sonic toothbrush

PHILIPS
sense and simplicity

Specifications

Additional features

- Exclusive patented Sonic technology: Dynamic fluid action helps clean between teeth
- Smartimer: Dentists recommend a two-minute brushing time
- 2 contoured brush head sizes: Standard for complete clean, Compact for precision
- 2 speeds: Choose the most comfortable speed for you
- Easy-Start: Gently increases power over 14 uses
- Programmable Quadpacer: 30-second interval timer for an even clean
- Deluxe recharge gauge: Green tracks battery charge, yellow for recharge
- Ergonomic tapered handle: Non-slip grip designed to improve handling

Technical specifications

- Voltage: 220-240

- Operating time (full to empty): 40
- Charging time: 24 hour(s)
- Battery rechargeable: NiCd
- Power consumption: 3 W
- Speed (top): Up to 31,000 movements per minute

Accessories

- HX7001: Standard Elite single pack brush head
- HX7002: Standard Elite double pack brush head
- HX7011: Mini Elite single pack brush head
- HX7012: Mini Elite double pack brush head

Weight and dimensions

- F-box dimensions: 17.2 x 9.5 x 23.5 in
- F-box weight: 680 lb
- A-box dimensions: 35.1 x 29.7 x 24.1 in
- A-box weight: 4.08 lb
- A-box number of pieces: 6 pcs

Highlights

Safe and gentle

Your Sonicare is safe to use on: braces (brush heads wear out sooner when used on braces), dental restorations (fillings, crowns, veneers) and periodontal pockets.

Healthier gums

Clinically proven to reduce gingivitis and help shrink periodontal pockets.



Issue date 2010-12-17

Version: 2.1.1

12 NC: 0000 000 00000

© 2010 Koninklijke Philips Electronics N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips
Electronics N.V. or their respective owners.

www.philips.com