



Philips
Digital Voice Tracer

512MB
33 h recording time



LFH0600

Notetaking at its best

Just remember

The Philips Digital Voice Tracer 600 is a digital recorder for those who think on their feet. Just start talking and the Digital Voice Tracer will do the rest. So why search for a pen when you can take accurate notes at the speed of sound?

Fast and simple

- Extra large display for an at-a-glance overview
- Instant one-touch recording - from power off to record mode
- Organize your files in a fast and easy way with four folders
- USB mass storage ensures maximum compatibility

Advanced

- Record in the popular MP3 format
- Rapidly transfer recordings and data via high-speed USB 2.0
- Voice-activated recording for hands-free note taking

Designed around you

- Lightweight and ergonomic design
- Extra large front speaker for crystal-clear voice playback
- The Hold function prevents unintentional on/off switching

PHILIPS

Specifications

Connectivity

- Headphone: 3.5 mm
- Microphone: Built-in microphone, Microphone socket
- USB: High-speed USB 2.0

Picture/Display

- Type: LCD
- Diagonal screen size: 3.5 cm / 1.5 inch

Storage Media

- Built-in memory (RAM): 512 MB
- Built-in memory type: NAND Flash
- Mass storage class compliant

Audio Recording

- Recording formats: MPEG1 layer 3 (MP3)
- Recording modes: HQ (MP3/mono), SP (MP3/mono), LP (MP3/mono)
- Bit rate: 64 kbps (HQ mode), 48 kbps (SP mode), 32 kbps (LP mode)
- Recording Time: 33.5 hours (LP mode), 22 hours (SP mode), 16.5 hours (HQ mode)
- Sample Rate: 22 kHz (HQ mode), 16 kHz (SP/LP mode)

Speaker

- Built-in speaker: 30 mm

- Speaker Type: Round dynamic speaker
- Output Power: 110 mW

Power

- Batteries included
- Battery type: AAA / LR03 Alkaline
- Battery lifetime: 30 hours (LP mode)

System Requirements

- PC OS: Windows 2000 / XP / Vista, Mac OS X, Linux
- USB: Free USB port

Green Specifications

- Lead-free soldered product

Accessories

- Included accessories: USB cable, batteries, user manual
- Optional accessories: Microphone/earphone combi, Plug-in Microphone LFH9171, Tie Clip Microphone LFH9173, Stereo Headphones LFH2236

Product dimensions

- Product dimensions (W x H x D):
4 x 11.8 x 1.8 cm
- Weight: 0.06 kg

Highlights

Voice-activated recording

Voice-activated recording is a convenience feature for hands-free recording. When enabled, recording will start automatically when you start speaking. Recordings stop automatically after you stop speaking and restart automatically when you start speaking again.

Extra large display

The large graphical display lets you view all important file information at a glance.

Instant one-touch recording

Press and hold the record button for a moment to switch the Digital Voice Tracer on and start recording instantly.

Easy file organization

The VoiceTracer digital recorder supports up to 99 recordings in each of four available folders to organize your files in a fast and easy way.

Lightweight design

Small and lightweight, you'll hardly notice it's with you. Pop it in your pocket or bag and take it anywhere. The durable casing and attractive, ergonomic design will make you want to carry it everywhere you go.

USB mass storage device

As a USB mass storage device, the Voice Tracer offers a convenient way to save, back up and move files. The Voice Tracer connects to a computer using a USB connection and does not require special software to be installed. The Voice Tracer is automatically displayed as a removable drive, allowing you to simply drag and drop files to and from the unit.



Issue date 2018-03-30

Version: 2.2.2

12 NC: 8735 060 00011
UPC: 0 37849 98349 2

© 2018 Koninklijke Philips N.V.
All Rights reserved.

Specifications are subject to change without notice.
Trademarks are the property of Koninklijke Philips N.V.
or their respective owners.

www.philips.com