

# PHILIPS

## Infrared lamp

300 W

For deep, zone treatment

Adjustable angle



PR3120/00



## Effective\*, drug-free pain relief

### Clinically proven technology\*

The InfraCare technology provides comforting warmth by penetrating deep into the skin. This increases local blood circulation & enhances the supply of oxygen & nutrients to help treat sore muscles and effectively recover from joint pain

#### **Effective\* pain relief**

- Pain relief and better mobility
- Helps the body to heal in a natural manner

#### **Proven health benefits**

- Clinically proven
- Classified medical device

#### **Zone treatment**

- Localized pain relief

# Highlights

## Effective\* pain relief




Our clinically proven technology provides relief from muscle and joint pains without the use of drugs. The InfraCare PR3120 stimulates local blood circulation and increases the local supply of oxygen for an effective zone treatment. This in turn provides a relaxing affect to the injured area and helps in an effective\* recovery from pain.

## Deep infrared treatment



InfraCare is designed to deliver the highest level of infrared waves reaching the deepest levels of the skin. Infrared warmth increases the local blood circulation, which speeds up the removal of metabolic end products and carbon dioxide from muscle tissue.

## Clinically proven



Trusted technology used by professionals. Clinically proven to objectively lead to functional joint improvement\*. InfraCare has proven to lead to significant pain reduction in 71% of rheumatoid arthritis patients and in 72% of low back pain patients\*.

## Classified medical device



A safe and effective way to treat your pain. InfraCare complies with IEC 60601-1 3rd edition for medical devices.

## Adjustable angle



The InfraCare 300W Lamp can be easily positioned from 0-40 degrees to relieve pain in a comfortable and an effective\* manner.

# Specifications

## Easy to use

Adjustability: 0-40 degrees (backwards)  
1 handle + 1 handgrip  
Mechanical timer

## Focused treatment

Treatment area: 40 x 30 cm  
Infrared halogen lamp: 300 W

## Logistic data

Country of origin: Hungary

## Medical appliance

Medical Device Directive: 2007/47/EC, MDD 93/42/EEC

## Safety

IEC certified: IEC 60601-1, 3rd edition  
30 minutes auto shut-off

## Technical specifications

Power: 300 W  
Voltage: 110 V or 220-230 V  
Frequency: 50-60 Hz  
Insulation: Class II (double isolation)  
Lifetime of lamp: 500 hour(s)  
Cord length: 250 cm

## Weight and dimensions

Product dimensions: 18,0 x 20,6 x 31,6 cm  
Product weight: 2 kg  
F-box dimensions: 22,3 x 24,9 x 33,4 cm

