

# PHILIPS

4 Self-adhesive  
electrodes

PulseRelief

For use with PulseRelief

4 electrodes

Patented hydrogel technology

Easy magnet connection

PR3827/00



## Muscle and Joint pain relief

### Self-adhesive electrodes

Philips self-adhesive electrodes have patented hydrogel technology for multiple application to the skin. There are 4 in a pack, and each are recommended for up to 20 uses. The electrodes are easily connected to the TENS Units by a magnet.

#### **Patented hydrogel technology**

- Patented hydrogel technology for multiple use

#### **Optimized for multiple use**

- Can be used for several treatment sessions, up to 20 times

#### **Magnet for easy connection**

- Easy magnet connection to connect TENS Units with electrodes

# Highlights

## Multiple use

Our electrodes use a high quality patented hydrogel technology for multiple application to the skin. The electrodes are flexible so they can be applied on irregular body contours.

## Reusable up to 20 times

Can be used for several treatment sessions (up to 20 times) on different body parts. The number of times the electrodes can be used is

dependent on skin type and usage. The electrodes can easily be stored in their original plastic package for next time use.

## Easy snap-on connection

The TENS units are easily connected to the self-adhesive electrodes by a magnet. First, apply the self-adhesive electrodes to your body and then just snap-on the TENS units via a magnet connection.

# Specifications

## What's in the box

Self-adhesive electrodes: 4 pcs

Self-adhesive electrodes size: 50 x 50 mm

Multistick<sup>®</sup> hydrogel

Plastic seal bag

## Composition of hydrogel

Glycerin, non-tallow: (Cas. No.: 56-81-5)

Water: (Cas. No.: 28062-44-4)

Poly(acrylate) Co-polymer: (Cas. No.: 7732-18-5)

