

PHILIPS

Steam generator
iron

PerfectCare 9000
Series

ActiveSense technology

Guaranteed no burns*

Light weight iron

Up to 750g steam boost



PSG9050/26



The only iron that knows what you are ironing Effortless ironing. Excellent results.

PerfectCare 9000 Series comes with the world's first fabric sensing technology called ActiveSense. Using a built-in camera and AI, it knows what you are ironing and adjusts the ideal temperature and steam amount for effortless results.

Effortless ironing with excellent results

- The only steam generator that knows what you're ironing
- Automatically adjusts settings
- Guaranteed no burns on all ironable garments

Easy crease removal with automatic powerful steam

- Intelligent automatic steam for faster and easier ironing
- Powerful steam for ultimate crease removal
- Light weight iron ideal for vertical steaming
- SteamGlide Elite Soleplate, ultimate gliding & durability

Extra convenient to use

- Large 1.8L detachable watertank for easy refilling
- Steam sound dampening technology
- Easy De-calc Plus
- Automatic shut off when iron is left unattended
- Safe carry lock for safe and easy transport

Highlights

Knows what you are ironing



ActiveSense technology recognizes fabrics using a built-in camera and artificial intelligence.

Automatically adjusts settings



It automatically adjusts temperature and steam amount to give excellent results for all fabrics.

Intelligent automatic steam



Our DynamiQ smart sensor knows exactly when and how the iron is moving, delivering powerful steam where you need it most – for easier ironing and faster results.

Ultra -powerful steam



Strong, continuous steam tackles even the thickest fabrics with ease. Watch stubborn creases melt away with an extra steam boost just where you need it.

No burns - guaranteed



Thanks to the fast cooling system, we guarantee this iron will never cause burns to any ironable fabric when it is left unattended. You can leave it resting on your clothes or ironing board.

Light weight iron



The iron is amazingly light and comfortable to handle, gliding easily and reducing wrist strain. Its minimal weight also makes it easy for vertical steaming of curtains and hanging clothes.

1.8L detachable water tank



1.8-liter transparent tank gives you up to 2 hours of continuous use. When the water tank is empty, you will be reminded by the indicator light and you can refill easily at any time under the tap through the large filling door.

SteamGlide Elite Soleplate



SteamGlide Elite is our best soleplate for ultimate gliding performance and maximum scratch resistance. Its stainless steel base is twice as hard as a regular aluminum base, and our patented 6-layer coating with advanced titanium layer glides effortlessly for the fastest results.

Easy De-calc Plus



Regular descaling protects your iron, extends its lifetime, and ensures the best steam performance. Our exclusive Easy De-Calc Plus system collects limescale continuously with an indicator light and audio signal to tell you when it needs emptying. Simply remove the plug and let the water and scale particles flow.

Specifications

Technology

ActiveSense technology

DynamiQ

Fast cooling soleplate

Steam sound dampening

Silent steam Technology

Fast crease removal

Power: 3100 W

Easy to use

Water tank capacity: 1800 ml

Detachable water tank

Hose length: 1.8 m

Power cord length: 1.8 m

Safety auto off

Heat up time: 2 min

Scale management

Descaling and cleaning: Easy De-calc Plus

Size and weight

Packaging dimensions (WxHxL): 33.4 x 35 x 57.4 cm

Product dimensions (WxHxL): 22.9 x 28.8 x 49.3 cm

Total weight with packaging: 8.8 kg

Weight of iron: 0.97 kg

Weight of iron + base: 5.5 kg

Green efficiency

Energy saving*: 22 %

General specifications

ActiveSense

AutoSteam with DynamiQ

Calc Clean reminder: Light & Sound

Cord storage: Compartment

Hose storage: Compartment

Max bar pressure: 9.0

Power cord material: Textile

Safety carry Lock

Soleplate: SteamGlide Elite

Steam boost up to (g): 750

Steam rate up to (g/min): 180

Ultra-light iron

Warranty: 2 + 1 years with registration

