

# PHILIPS

Adjustable crank set

ReActiv

PTE4010CA



## Upgrades bikes to fit all users

Best for folks with limited mobility

The Adjustable crank set makes exercises in smaller range of motion possible on bikes. This accessory gives greater accessibility to users with needs to improve range of motion and helps reduce fear of pain.

**Highly adjustable for different range of motion**

- This accessory enables closed kinetic chain exercises.

# Highlights

## Adjustable crank set

The crank is incrementally adjustable to adapt to a user's limited range of motion, whether the user has just resumed physical activities or has made good progress towards recovery. Without being far along in the recovery

process, users can start exercising with just 15 degrees range of motion, thereby reducing the fear of pain. Users can gradually increase range of motion when they feel comfortable to complete a 360-degree pedaling range.

# Specifications

## Accessory set

set: 1

## Balance

-: Weighted

## Product material

-: ABS plastic

## Thread size

inches: Standard 9/16" thread

## FK-Physical therapy solutions

Types: ReActiv

