

PHILIPS

Recumbent bike

ReCare



PTE7000MR



Strengthen knees and ankles without the impact

With real time feedback

The 7.0 R Recumbent bike sets the standard for providing both comfort and effective exercises. The rotating seat allows users to easily get on and off the bike. Patients can be positioned properly for exercises.

Accessibility and comfort combined

- Adjust crank length to fit range of motion
- The seat rotates 360 degrees and adjusts forward and back

User friendly

- The display is easy-to-read for metrics
- The seat and pedal adjustments are numbered

Innovative programs

- Console measures strength input on left and right sides
- Displays workout data

Highlights

Console



Large, bright LED windows are easy for patients and clinicians to view important workout data at a glance.

Crank adjustment

To help patients with limited mobility or injuries, the pedals can be moved along the crank to accommodate range of motion as small as 15 degrees. Smaller degrees render smaller circular lower body motions, so that patients can complete a full 360 degree stroke during rehabilitation. Cranks are indexed for keeping track of how far patients have improved, as well as for precise repetition.

Indexed adjustments



Patients can document their bike set up using the indexed seat distance and crank adjustments. With documented numbers, you can easily reposition these features suitable for your next exercise or rehabilitation session.

Rotatable seat and adjustments

The padded seat features contoured curves for the best stability and comfort. This helps people with limited strength and those who need assistance to position themselves better while seated. The seat moves forward and back to accommodate different heights, as well as reclines for various hip angles.

Symmetry program

For asymmetrical patients who need to work on strengthening a weaker leg, the symmetry program allows clinicians to keep record and measure improvement as their patients progress.

Console

Patients and clinicians can quickly see essential, real-time performance feedback. Real-time streaming of patients' data is also available for clinicians to download. Windows display time, rotations per minute, watts, calories, METs, heart rate, and power.

Specifications

Compliance

Certifications: CAN/CSA-C22.2 No. 60601-1:14, ANSI/AAMI, ES60601-1:2005+A2 (R2012) +A1, IEC 60601-1-2:2014, EN 60601-1-2:2015, IEC 60601-1:2015+A1:2012, EN 60601-1:2006, +A1:2013, +A12:2014, EN ISO 20957:2013, MDD 93/42/EEC Class Im., NB: 0123

Max user weight

lbs: 440

Net weight

lbs: 180

kg: 82

Overall dimensions

inches: 55.1x32.6x46.8

cm: 140x83x119

Power AC

Volts: AC 100-240V

Drive train

Heavy-duty: 8 groove poly-V belt, with spring-loaded idlers

CSAFE

Yes: -

Hand pulse

Yes: -

Programs

include: Custom, Hill, HR, Interval, Manual, Plateau, Symmetry, VO2, Isokinetic

Resistance

level: 1 to 50

type: Magnetic resistance system

Start resistance

watts: 5

Work range

watts: 5 to 750

