

# PHILIPS

Hand / wrist  
stabilization set

ReCare



PTE7030MA



## Added stability for users

Help to get stable and correct hand position

This accessory allows patients with weaker hand or wrist to experience upper body conditioning with a firm grip. With reduction in pressure points, neurological patients with little hand strength can enjoy the benefits of ReCare steppers.

### **More accessibility for neurological patients**

- Glove-like fit secures hand to handle

# Highlights

## Hand wrap



Padded for comfort, the hand wrap secures hands with hook and loop fastener. It fits like a glove and helps to keep hands in proper positioning.

# Specifications

## Accessory set

set: 1

## Adjustment

-: Hook and loop straps

## Product material

-: Neoprene foam rubber

